

ICP-MS Automation Script

Code Author: Mohamed Harhash, PhD student in Environmental Analytics.

Primary Email: m.harhash@daad-alumni.de

Alternative Email: Harhash007@yahoo.com

Acknowledgment

This code was originally developed by Mohamed Harhash, under the supervision of Dr. Lars Duester, during his Ph.D. at the Federal Institute of Hydrology, Koblenz, Germany.

If you build upon or further develop this code, please acknowledge or cite the original work. For any inquiries, you can contact Mohamed Harhash at the provided emails.

Citation

Harhash, M. (Supervised by Dr. Lars Duester, Federal Institute of Hydrology, Germany). (2024). *Automation solution for ICP-MS* (Version 1.0; Python code). GitHub. Available at <https://github.com/Harhash2024/ICP-MS-AUTOMATION.git>

Objective:

The main objective was to automate an inductively coupled plasma mass spectrometer (ICP-QMS 7700 / ICP-QQQ- MS 8900) for the monitoring of metal(loids) in river water in a single measurement run with a high time resolution of one hour. To achieve this objective, a self-constructed sample collector/autosampler was constructed to collect samples continuously, 24 hours a day, and an automated ICP-MS instrument was programmed to automatically initiate and perform measurement sequences at scheduled times.

Code capabilities :


The code is developed to automatically initiate and manipulate the ICP Go software to operate the Agilent ICP-MS 7700 and ICP-QQ-MS 8900 instruments or any other Agilent instrument compatible with the ICP Go software. It carries out measurement sequences and automates the process : 1) logging in to the ICP Go software, 2) turning on the plasma, 3) waiting for warm-up time, 4) selecting and creating batches, 5) performing batch tuning, 6) starting batches and measurements, 7) monitoring errors, 8) sending emails about instrument status, and after measurements, 9) turning off the plasma, 10) logging out and closing the software.

Content

Before you start:.....	3
Instructions for using ICP Go software	3
Python automation code:.....	5
Steps to run the code	5
Code Explanation with Snippet	8
1. Open the ICP Go software	8
2. Logging into the ICP Go software.	9
3. Turning on the plasma.....	10
4. Creating a batch.....	11
5. Adding the batch to the queue.....	13
6. Warming up time.....	14
7. Run the queue and start tuning and the batch.....	15
8. Waiting tuning of the batch till finish	16
9. Check the error and send an email about the instrument status	17
10. Waiting time for the measurement	18
11. Log out and close ICP Go software.....	19

Before you start:

- Make sure that the autosampler software is always opened and minimized.



- ICP Go software is an interface to operate and control Agilent ICP-MS machines, and you have to run ICP Go software in the background when running automation code.

Instructions for operating ICP Go software

1. Double-Click to Initiate ICP Go Software:

- Locate the ICP Go software icon  on your Desktop.
- Double-click the icon to start the program.

2. Running in the Background:

- Once the ICP Go initiated, the software will continue to run in the background, and It will remain active until you manually close it.

Make sure that the Masshunter service and ICP Go software is running in background





3. Accessing the ICP Go Software Interface:

- Double-click the ICP Go icon on your computer. This will automatically open your web browser to the correct URL (<http://xxx.xxx.xxx.xxx:5000>) to access the ICP Go interface.
- If needed, you can also copy this **URL** from ICP Go and paste it into any web browser's address bar to open the software manually.

- The ICP Go user interface will open in the browser window, allowing you to control the ICP-MS machine. Save this URL from ICP Go and replace it in the ICP-MS automation code.

4. Operating the ICP-MS machine:

- Use the web-based interface ICP Go software to perform tasks such as turning plasma on, creating batches, auto-tuning, starting analyses, monitoring instrument status, and adjusting settings. It has all control options which are provided by MassHunter software.



Python automation code:

Steps to run the code

1. Install python 3.x

To install Python through Anaconda:

- Go to the official anaconda website and download it.
- Choose the version suitable for your operating system, and anaconda includes Python 3 by default.
- Windows: Install anaconda by double click `anaconda.exe` file and follow the instructions.
- After installation, you can verify installation by opening Anaconda **prompt** and typing this command `conda --version` to know the conda version.
- To keep python projects organized, create separate environment for each project.

```
Conda create -n myenv python=3.9    ##### you can replace myenv with the name of your
environment and your python version.
```

- Activate the environment
`conda activate myenv`
- Launch Python or Jupyter
 - You can do this through the Anaconda interface.
 - To start a Jupyter Notebook, click the Launch button next to Jupyter Notebook.

2. Download the code:

- a) You can download** the code directory from this GitHub repository.

<https://github.com/Harhash2024/ICP-MS-AUTOMATION.git>

OR

- b) Clone the repository:**

Install Git:

Download Git from the official website: <https://git-scm.com/downloads>.

Run the installer and follow the prompts to complete the installation.

Restart your **command prompt** to ensure the new git command is recognized.

Command prompt:

```
git clone https://github.com/Harhash2024/ICP-MS-AUTOMATION.git
cd ICP-MS-AUTOMATION
```

1) Install Python dependencies:

Install the necessary Python packages using the `requirements.txt` file:

```
pip install -r requirements.txt
```

3. Set up ChromeDriver: Download ChromeDriver from

<https://developer.chrome.com/docs/chromedriver>

Ensure that the ChromeDriver is either in your system's `PATH` or provide the full path in the script:

```
python
driver =
webdriver.Chrome(executable_path='C:/path/to/chromedriver.exe')
```

4. Configure Email Notifications:

Update the script with your email credentials to enable the error notification feature:

```
python
gmail_email = "your-email@gmail.com"
gmail_password = "your-app-password"
recipient_email = "recipient-email@example.com"
```

5. Ensure ICP Go Software Accessibility:

Ensure that the ICP Go software is accessible via the web interface at `http://xxx.xxx.xxx.xxx:5000`. ***Update this in the code as well.***

To run the script, navigate to the directory containing it and execute the following steps:

1. Place your Python script (`code_name.ipynb`) in a preferred folder.
2. To schedule and run the code, create a `.bat` file (a `.bat` file is used to execute commands in Windows).
3. In the notebook file, copy and paste the following command lines, modifying the paths to match the locations of your folders.
 - o Command lines:

```
@echo off
```

```
:: Path to the Jupyter executable
"C:\path\to\jupyter.exe" nbconvert --to notebook --execute
"C:\path\to\your_notebook.ipynb" --output
"C:\path\to\output_notebook.ipynb"
```

- Name this notebook file as you prefer e.g., **ICP-MS_Automation.bat**

4. Schedule the .bat File in Task Scheduler

- Press **Windows Key + R** on your keyboard to open the **Run** dialog box.
- Type **taskschd.msc** and press **Enter**. This will open the **Task Scheduler**.
 - a. Open **Task Scheduler** and create a new task.
 - b. In the **Actions** tab, select **Start a program**.
 - c. Browse to the .bat file you created (e.g., **ICP-MS_Automation.bat**).
 - d. Set the **Triggers** and other settings as needed, then save the task.


With this setup, your .bat file will directly execute your Jupyter Notebook file each time it's run.

Script Workflow:

1. **Login:** The script logs into ICPMS Go software using stored credentials in the code.
2. **Plasma Control:** The plasma is turned on.
3. **Batch Creation:** The script creates and saves a new batch using the appropriate template based on the current date and time.
4. **Warm-up:** It waits for the plasma to warm up.
5. **Batch Queueing:** Adds the created batch to the queue for measurement.
6. **Tuning:** Performs any necessary tuning for the measurement.
7. **Measurement:** The script waits for the measurement to complete.
8. **Error Handling:** If any errors are encountered, an email notification is sent, and the error is logged.
9. **Turn off the plasma**
10. **Logout and Shutdown:** Once the batch is complete, the script logs out of ICPMS Go and safely closes the software.

Code Explanation with Snippet

1. Open the ICP Go software




Code snippet:

```
#1- Opening ICPMS Go software

driver = webdriver.Chrome()
driver.get("xxx.xxx.xxx.xx:5000")
### Check this IP address if it is the same for your ICP-MS. it is provided
with the software
driver.maximize_window()
```

2. Logging into the ICP Go software.

- You can modify **user name** and **password** according to your system



Code snippet:

```
#2-Logging into ICPMS Go
wait = WebDriverWait(driver, 20)
username_input = wait.until(EC.element_to_be_clickable((By.NAME,
"username")))
username_input.click()
username_input.send_keys("your_username") ##### replace it with your own user
name
password_input = driver.find_element(By.NAME, "password")
password_input.send_keys("Mondel") ##### replace it with your own password
login_button = driver.find_element(By.ID, "submit")
login_button.click()
## waiting 7 second to laod the page
time.sleep(7)
...
```


3. Turning on the plasma

The screenshot shows a 'Batch Management' interface with a grid of 12 batch entries. Each entry includes a delete icon. The first row contains four entries: 'Wednesday morning_18_09_20...', 'Tuesday evening_17_09_2024 (36/36)', and two instances of 'Tuesday morning_17_09_2024 (35/35)'. The second row contains four entries: 'Monday evening_16_09_2024 (35/35)', 'Monday morning_16_09_2024 (35/31)', 'Sunday morning_15_09_2024 (35/35)', and 'Saturday evening_14_09_2024 (34/31)'. The third row contains four entries: 'Saturday morning_14_09_2024 (35/31)', 'Friday evening_13_09_2024 (35/35)', 'Friday morning_13_09_2024 (35/35)', and 'Thursday evening_12_09_2024 (36/31)'. The fourth row contains two entries: 'Friday evening_13_09_2024 (35/35)' and 'Thursday evening_12_09_2024 (36/31)'. A blue box highlights the 'Plasma button' in the top right corner, which has an arrow pointing to a water droplet icon.

Code snippet:

```
#3-Starting the plasma
driver.find_element(By.CSS_SELECTOR, "#plasma .pe-ripple__mask").click()
time.sleep(2)
driver.find_element(By.CSS_SELECTOR, "#yes .pe-ripple__mask").click()
```

4. Creating a batch



Code snippet :

```
#4-Creating and managing new batches based on templates
driver.find_element(By.CSS_SELECTOR, "svg").click()
driver.find_element(By.ID, "template").click()


# determine the current date, time and day of the week
current_day_index = datetime.datetime.now().weekday()
current_time = datetime.datetime.now().time()
current_time_period = "morning" if current_time <
datetime.datetime.strptime("12:00:00", "%H:%M:%S").time() else "evening"
day_names = ["Monday", "Tuesday", "Wednesday", "Thursday", "Friday",
"Saturday", "Sunday"]
selection = f'{day_names[current_day_index]}_{current_time_period}'

# click dropdown list for selecting proper batch for the current day and time
time.sleep(4)
dropdown_element = driver.find_element(By.ID, "template")
dropdown_element.click()
dropdown = Select(dropdown_element)
dropdown.select_by_visible_text(selection)

# Naming the created new batch as day name 'evening' or 'morning' and date
timestamp = datetime.datetime.now().strftime("%d_%m_%Y")
element_name = f'{selection}_{timestamp}'
print(element_name)
driver.find_element(By.NAME, "batchname").click()
time.sleep(5)
driver.find_element(By.NAME, "batchname").send_keys(element_name)
time.sleep(10)

# Saving the batch
driver.find_element(By.CSS_SELECTOR, "#savebatch .pe-ripple__mask").click()
time.sleep(10)
```

- Open the list of batch templates and select the appropriate batch.



Dynamic creation and management of batches


- As the automation is designed for daily monitoring, the ICP-MS machine will run twice a day, morning and evening, to measure multiple elements in river water. Therefore, two batch templates are created for each day of the week.
- Each batch is configured for the appropriate location of the sample vials and the two blanks based on the corresponding day.
- The evening batch is configured to measure samples in Rack 1 and Rack 2, which contain samples collected between 7:00 and 18:00.

- The morning batch is configured to measure samples in Rack 3 and Rack 4, which contain samples collected between 19:00 and 6:00 on the previous day.
- For clarification, a sample named '7:00' represents the mixed sample collected from 7:00 to 8:00.

Naming System of created batch.

The code will first determine the current date, time, and day of the week (e.g., Monday, Tuesday...). It will then identify whether it is morning or evening, after 12:00 is considered evening, and before 12:00 is considered morning. Based on this information, the code will select the appropriate batch, such as "Monday evening" from the batch list, and name it as 'Monday evening _ date of measurement'.


5. Adding the batch to the queue




Code snippet :

```
#Add the created batch to queue

time.sleep(25)
xpath_selector = f'//*[@id="queuebutton_{element_name}"]/div'
element = driver.find_element(By.XPATH, xpath_selector)
element.click()
```



6. Warming up time



Code snippet :

```
#6- waiting warming time
time.sleep(1900)
```

7. Run the queue and start tuning and the batch

The screenshots show the 'Batch Management' screen of the MondeICP Go software. The interface includes a sidebar with icons for Home, Settings, and Logout. The main area displays a grid of batch entries:

- Wednesday evening_18_09_2024 (0/ X)**: Status bar shows 'Wednesday evening 2024-09-18T17:00:44'.
- Wednesday morning_18_09_20... (3! X)**: Status bar shows 'Wednesday morning 2024-09-18T05:00:44'.
- Tuesday evening_17_09_2024 (36/36 X)**: Status bar shows 'Tuesday evening 2024-09-17T17:00:44'.
- Tuesday morning_17_09_2024 (35/3 X)**: Status bar shows 'Tuesday morning 2024-09-17T05:00:44'.
- Monday evening_16_09_2024 (35/35 X)**: Status bar shows 'Monday evening 2024-09-16T17:00:42'.
- Monday morning_16_09_2024 (35/3! X)**: Status bar shows 'Monday morning 2024-09-16T05:00:44'.
- Sunday morning_15_09_2024 (35/35 X)**: Status bar shows 'Sunday morning 2024-09-15T05:00:44'.
- Saturday evening_14_09_2024 (34/3- X)**: Status bar shows 'Saturday evening 2024-09-14T17:00:44'.
- Saturday morning_14_09_2024 (35/3! X)**: Status bar shows 'Saturday morning 2024-09-14T05:00:44'.
- Friday evening_13_09_2024 (35/35) X**: Status bar shows 'Friday evening 2024-09-13T17:00:44'.
- Friday morning_13_09_2024 (35/35) X**: Status bar shows 'Friday morning 2024-09-13T05:00:44'.

In the bottom screenshot, a tooltip 'Pause the current running batch (Disabled)' is visible over the 'Tuesday evening_17_09_2024' batch card. A progress bar at the bottom indicates 'Wednesday evening_18_09_2024 progress 0/36 (1 batches in queue)'.


Code snippet :

```
# 7- Run the queue and start tuning and the batch
```

```
run_button_selector = "#run .pe-ripple__mask"
queue_button_selector = "#queue .pe-ripple__mask"

click_with_retry(driver, run_button_selector)
time.sleep(20)
click_with_retry(driver, queue_button_selector)
```

8. Waiting tuning of the batch till finish



Code snippet :

```
#8- Waiting tuning of the batch till finish
..
time.sleep(180)
..
```

9. Check the error and send an email about the instrument status

The screenshot shows the 'Instrument' tab of the MondelCP Go software. On the left, there's a sidebar with icons for Home, Favorites, and a gear (Settings). The main area has two sections: 'Hardware Configuration' (showing an Autosampler and an Instrument Test Report) and 'Hardware Information'. Under 'Hardware Information', there are four tables: 'Mainframe' (IP 192.168.1.128, FW 5.02.03 001.03, Model G3665A, Serial SG23233242, Software D.1.2.708.6); 'Sample Introduction' (Autosampler CETAC ASX520, Sample Introduction ISIS 3, Nebulizer MicroMist); 'Plasma' (Ignition Mode Aqueous Solution); and 'Ion Lenses' (Model x-Lens). Below these is an 'Error Log' section with a table:

ID	Severity	Time Stamp	Message

The top screenshot shows the Microsoft Outlook inbox. It includes a ribbon with File, Home, Send / Receive, Folder, View, Help, and a search bar. The left sidebar shows folders like Favorites, Posteingang, and Sent Items. The inbox lists two items: an unread message from 'Sonnet, Susanne' and a new message from 'harhash008@gmail.com'.

The bottom screenshot shows a search results page for 'error'. The search bar at the top contains 'error'. The inbox lists several messages from 'harhash008@gmail.com' with subject lines like 'Status_ICPMS_Monde_MH' and 'Status_ICPMS_Monde_MH'. The right side of the screen shows the message preview pane with the first message's content:

Error Message: Error Log
ID Severity Time Stamp Message

1 ERROR 2024-05-22 05:09:20 [Error] Disconnected from ICP-MS MassHunter Service provider:There was no endpoint listening at net.pipe://localhost/ICPMS_MASSHUNTER_MHAnalyzerInterface/IcpMHServiceProvider that could accept the message. This is often caused by an incorrect address or SOAP action. See InnerException, if present, for more details.

Code snippet :

```
#9- check the error and send an email about the instrument status

try:
    instrument_button = WebDriverWait(driver, 10).until(
        EC.presence_of_element_located((By.CSS_SELECTOR, "#instrument .pe-
ripple__mask"))
    )
    instrument_button.click()
except Exception as e:
    print(f"Failed to click the instrument button: {e}")

error_message_selector = "#errorLog > div"
try:
    error_message_element = WebDriverWait(driver, 10).until(
        EC.presence_of_element_located((By.CSS_SELECTOR,
error_message_selector))
    )
    error_message_text = error_message_element.text
    if error_message_text:
        gmail_email = "your_email@gmail.com" ##### your E-mail
        gmail_password = "xxx xxx xxx xxx" ## access code of your E-mail
        obtained from you email setting ## your-app-password
        recipient_email = "recipient-email@example.com" ### here you can use
your E-mail
        subject = "Error Alert_ICPMS_Monde" #####
        body = f"Error Message: {error_message_text}"

        msg = MIMEText(body)
        msg["Subject"] = subject
        msg["From"] = gmail_email
        msg["To"] = recipient_email

        server = smtplib.SMTP("smtp.gmail.com", 587)
        server.starttls()
        server.login(gmail_email, gmail_password)
        server.sendmail(gmail_email, recipient_email, msg.as_string())
        server.quit()
except Exception as e:
    print(f"An error occurred: {e}")

print("Email sent successfully.")
```

10. Waiting time for the measurement


- You can customize your needs

Code snippet :

```
# 10- Waiting time for the measurement (around 4 hours) and it is more than
our measurement time to ensure enough time after measurement ####

driver.find_element(By.CSS_SELECTOR, "#queue .pe-ripple__mask").click()
time.sleep(16000)
```

11. Log out and close ICP Go software.



Code snippet :

```
# 11- Sing out and close ICP Go software

try:
    logout_button = WebDriverWait(driver, 10).until(
        EC.element_to_be_clickable((By.CSS_SELECTOR, "#login .pe-ripple__mask"))
    )
    logout_button.click()
except TimeoutException:
    print("Logout button not clickable. Continuing without clicking.")
try:
    ok_button = WebDriverWait(driver, 10).until(
        EC.element_to_be_clickable((By.CSS_SELECTOR, "#ok .pe-ripple__mask"))
    )
    ok_button.click()
except TimeoutException:
    print("Ok button not clickable. Continuing without clicking.")

time.sleep(5)
driver.quit()
```

Alternative Methods to schedule the code initiation

Use a Robust Scheduling System:

- Use a robust scheduling system, such as Windows Task Scheduler, as it is managed by the operating system and can automatically restart failed tasks.
- You can set up triggers like 'on system start-up' or 'on user logon' to ensure the task restarts if it fails. Another advantage is that you don't need to initiate the code or Anaconda each time the computer restarts.

1. Using windows task scheduler:

- The Windows Task Scheduler allows you to schedule scripts, programs, or tasks to run at specific times or intervals.
- **Steps:**
 1. Open windows task scheduler.
 2. Click on "Create Basic Task" or "Create Task."
 3. Name the task and set a trigger (e.g., daily at a certain time).
 4. In the "Action" tab, select "Start a program" and browse to the script or executable you want to run.
 5. Set any conditions or settings as required and save the task.
- This method is reliable for scheduling on Windows machines and doesn't require changes to the code.

However, there are also other methods by which the code can be scheduled.

This includes:

2. Embedded Scheduling in the Code:

- You can include scheduling logic directly within the code itself using specific libraries or modules.

For example, in Python, you can use the schedule or time libraries to set specific times for the code to run.

3. Using Python Scheduler Libraries:

- Libraries such as AP Scheduler (Advanced Python Scheduler) offer flexible scheduling.